



DEVOLPE
MEDICAL

VACUUM SPLINTS

LEG



ARM



ANKLE



HEAD



REF

MANUFACTURER: DEVOLPE s.r.o.

REGISTERED TRADE MARK: DEVOLPE

REGISTERED ADDRESS: Nam. 14. října 1307/2, Prague 5, 150 00, Czech Republic

T/F: +420 605 700 078

Website: www.devolpe.com

Email: info@devolpe.com



DEVICE DESCRIPTION

A vacuum splint is used as a temporary splint in emergency situations. The thousands of polystyrene beads inside the splint work as an orthopedic cast by molding around the injured limb. The vacuum splint is ideal for stabilizing luxations, subluxations, and fractures. X-ray, CT, and MRI scans can be performed with the splint in place.

Material: material combination with PU, PVC, EVA and new plastic materials.

TECHNICAL CHARACTERISTICS

Product	Type	Models	Catalog code	Length (cm)	Width (cm)	Height (cm)	Weight (kg)
Vacuum Arm Splint	Adult	Solo-chamber	VS/03/AVSSDP	67	50/30	3	0.9
Vacuum Arm Splint	Adult	Poli-chamber	VS/04/AVSPDP	67	50/30	3	0.9
Vacuum Leg Splint	Adult	Solo-chamber	VS/01/LVSDP	110	70/40	3	1.8
Vacuum Leg Splint	Adult	Poli-chamber	VS/02/LVSPDP	110	70/40	3	1.8
Vacuum Ankle Splint	Adult	Solo-chamber	VS/05/ANVSSDP	57	47	3	0.9
Vacuum Ankle Splint	Adult	Poli-chamber	VS/06/ANVSPDP	57	47	3	0.9
Vacuum Head Immobilizer	Adult	Poli-chamber	VS/07/HVSDP	31	38	8	0.7

Made in accordance with: ISO 13485:2016 standard

Accessories



Manual Pump
For vacuum products



Vacuum splints
bag